

Supporting Information

Towards an automated approach for rapid separation of actinides using a liquid handling system

*Nicholas P. Richard**, Sean, R. Scott, Christian Berry, Matthew RisenHuber, William Munley,
Karen Noyes, Matt Douglas, Lori Metz

¹Pacific Northwest National Laboratory, Richland, WA, 99352 USA

Table of Contents

Figure S1. Image of the 3D printed column piece positioned in the 3D printed column holder with the 3D printed positive pressure piece attached to a compressed air line.....	2
Figure S2. Hamilton setup for microchemistry.....	3




Figure S1. Image of the 3D printed column piece positioned in the 3D printed column holder with the 3D printed positive pressure piece attached to a compressed air line.




Figure 2. Hamilton setup for microchemistry showing custom-designed and purchased components. The sample wells were custom-machined specifically for this project, while the columns, collection wells, column rack, and positive air pressure system were 3D-printed. The heat shaker, pipette tips and tip rack, waste container, chemical reservoirs, and collection well are commercial components purchased with the Hamilton system.