

Supplementary Information (SI) for Journal of Materials Chemistry B.

This journal is © The Royal Society of Chemistry 2025

# Highly Stable Nanosized Lucigenin@PS-COOH: A Robust Chemiluminescence System for Sensitive cTnT Analysis

Fugang Li,<sup>a</sup> Qiang Li,<sup>\*c</sup> Feng Chen,<sup>a</sup> Chenxu Yan,<sup>b</sup> Yongxi Zhao,<sup>a</sup> and Zhiqian Guo<sup>\*b</sup>

---

[a] F. Li, Prof. Dr. F. Chen, Prof. Dr. Y. Zhao


Institute of Analytical Chemistry and Instrument for Life Science,  
The Key Laboratory of Biomedical Information Engineering of Ministry of Education,  
School of Life Science and Technology,  
Xi'an Jiaotong University, Xi'an, Shaanxi, 710049, PR China  
E-mail: chenxuyan@ecust.edu.cn;

[b] Prof. Dr. Z. Guo, Prof. Dr. C. Yan

Key Laboratory for Advanced Materials and Joint International Research Laboratory of Precision Chemistry and Molecular Engineering,  
Shanghai Key Laboratory of Functional Materials Chemistry,  
Feringa Nobel Prize Scientist Joint Research Center, Institute of Fine Chemicals,  
Frontiers Science Center for Materiobiology and Dynamic Chemistry,  
School of Chemistry and Molecular Engineering  
East China University of Science and Technology, Shanghai 200237, China  
E-mail: guozq@ecust.edu.cn;

[c] Dr. Q. Li

National Key Laboratory of Intelligent Tracking and Forecasting for Infectious Diseases,  
National Institute of Parasitic Diseases at Chinese Center for Disease Control and Prevention,  
Chinese Center for Tropical Diseases Research, National Research Center for Tropical Diseases,  
Key Laboratory of Parasite and Vector Biology, National Health Commission,  
WHO Collaborating Center for Tropical Diseases, National Center for International Research on Tropical Diseases,  
Ministry of Science and Technology, Shanghai 20025, China  
E-mail: liqiang@nipd.chinacdc.cn.



**Fig S1.** Standard curve of chemiluminescence (CL) intensity vs. the concentration of the cTnT calibrated substance using AE.



**Fig S2.** UPPER CycloneSun chemiluminescence automatic immunoanalyzer.

**Table S1** Chemiluminescence (CL) intensities of lucigenin@PS-COOH and acridinium ester in response to increasing test values of cTnT calibrated substance.

	cTnT calibrated substance (ng/mL)	CL intensity of lucigenin@PS-COOH (a.u.)	CL intensity of acridinium ester (a.u.)
1	0	1106	813
2	0.046	11992	7864
3	0.454	151274	191368
4	1.523	596965	387524
5	4.269	1817255	2636882
6	10.315	4772407	3976754

**Table S2** Repeated measurements of the cTnT quality control substance used to determined assay sensitivity ( 0.32 ng/L) and coefficient of variation (CV, 7.77%).

	Test (ng/mL)
1	0.28
2	0.31
3	0.29
4	0.28
5	0.33
6	0.32
7	0.31
8	0.35
9	0.28
10	0.30
Average	0.305

**Table S3** Comparison of detection kit and this work

	detection limit	calibration range, R <sup>2</sup>	CV
detection kit	≤ 0.3 ng/mL	≥ 0.9900	≤ 10%
This work	0.25 ng/mL	0.9988	7.77%

**Table S4** Comparison of Longhua Hospital results and this work for 50 clinical serum samples. The difference for each sample was within 15%, and the overall difference did not exceed 15%.

Sample	Longhua (ng/mL)	This work (ng/mL)	Comparison
1	7.985	7.707	103.60%
2	3.540	3.264	108.47%
3	6.547	6.239	104.94%
4	5.018	4.707	106.61%
5	6.321	6.010	105.17%
6	5.234	4.921	106.36%
7	7.540	7.250	104.00%
8	5.497	5.183	106.06%
9	4.597	4.292	107.12%
10	2.750	2.512	109.49%
11	4.026	3.733	107.84%
12	7.005	6.704	104.49%
13	1.057	0.971	108.90%
14	8.295	8.027	103.34%
15	1.689	1.531	110.29%
16	2.526	2.301	109.75%
17	4.506	4.202	107.23%
18	4.521	4.217	107.21%
19	2.451	2.231	109.84%
20	4.698	4.391	106.99%

21	1.987	1.803	110.23%
22	3.009	2.756	109.17%
23	3.210	2.948	108.91%
24	0.889	0.826	107.66%
25	4.351	4.050	107.42%
26	1.240	1.131	109.67%
27	8.015	7.738	103.58%
28	0.705	0.670	105.28%
29	1.058	0.971	108.90%
30	0.594	0.577	102.92%
31	0.720	0.682	105.53%
32	0.654	0.627	104.31%
33	0.447	0.457	97.84%
34	0.689	0.656	104.99%
35	0.508	0.506	100.31%
36	0.298	0.338	88.08%
37	0.254	0.304	83.50%
38	1.241	1.132	109.67%
39	0.702	0.667	105.22%
40	0.825	0.771	106.99%
41	0.775	0.729	106.36%
42	0.651	0.624	104.25%
43	0.334	0.367	91.11%
44	0.497	0.497	99.91%
45	0.608	0.589	103.27%
46	2.999	2.747	109.18%
47	1.442	1.310	110.10%
48	0.278	0.323	86.13%
49	0.560	0.549	101.99%
50	0.560	0.549	101.99%