

## *Supporting Information*

# **Electrochemical pinacol coupling of aromatic carbonyl compounds in a [BMIM][BF<sub>4</sub>]/H<sub>2</sub>O mixture**

Hannah Kronenwetter, Jakub Husek, Brian Etz, Aaron Jones, and Renuka Manchanayakage\*

Department of Chemistry, Susquehanna University, 514 University Avenue, Selinsgrove, PA 17870

manchanayakage@susqu.edu

### **Contents**

General information and the procedure for pinacol coupling reaction .....	2
Cyclic voltammogram of 80% [BMIM][BF <sub>4</sub> ]/H <sub>2</sub> O .....	2
NMR spectra of 1,2-diphenylethane-1,2-diol [4a] .....	3
NMR spectra of 1,2-di(4-methylphenyl)ethane-1,2-diol [4b] .....	5
NMR spectra of 1,2-di(4-methoxyphenyl)ethane-1,2-diol [4c] .....	7
NMR spectra of 1,2-di(4-chlorophenyl)ethane-1,2-diol [4d] .....	9
NMR spectra of 1,2-di(4-bromophenyl)ethane-1,2-diol [4e] .....	11
NMR spectra of 2,3-diphenylbutane-2,3-diol [6a] .....	13
NMR spectra of 2,3-di(4-methylphenyl)butane-2,3-diol [6b] .....	15
NMR spectra of 2,3-di(4-methoxyphenyl)butane-2,3-diol [6c] .....	17
NMR spectra of 2,3-di(4-bromophenyl)butane-2,3-diol [6d] .....	19
NMR spectra of 2,3-di(2-pyridinyl)butane-2,3-diol [6e] .....	21
NMR spectra of 3,4-diphenylhexane-3,4-diol [6f] .....	23
NMR spectra of 2-methyl-1-phenylpropan-1-ol [7a] .....	25
NMR spectra of diphenylmethanol [7b] .....	27

## General Information

All reactions were conducted under nitrogen atmosphere. All reagents were obtained from commercial suppliers and used without further purification. Ionic liquids, [BMIM][BF<sub>4</sub>] and [BMIM][PF<sub>6</sub>] were synthesized according to the literature.<sup>1</sup> NMR experiments were conducted in chloroform-*d*.

## General procedure for the electrochemical pinacol coupling reaction

The carbonyl compound (0.35 mmol) and a mixture of 80% [BMIM][BF<sub>4</sub>]/H<sub>2</sub>O (4.0 mL of [BMIM][BF<sub>4</sub>] and 1.0 mL of H<sub>2</sub>O) were placed in an electrochemical cell fitted with a tin foil anode (1.0 cm<sup>2</sup>), a platinum plate cathode (1.0 cm<sup>2</sup>) and a Ag/AgCl reference electrode. The mixture was stirred and degassed by bubbling nitrogen for 30 minutes. A controlled potential of 2.0 V was applied under nitrogen atmosphere for 5 hours. The resultant solution was extracted with diethyl ether (2 × 5 mL). The ether extract was dried with anhydrous sodium sulfate, filtered and concentrated in vacuo. When necessary, the residue was purified via silica gel flash column chromatography (10% ethyl acetate in hexane eluent) to afford the final product.

## References

1. Dupont, J.; Consorti, C.; Suarez, P.; Souza, R. *Org. Syn. Coll. Vol.* **2004**, 79, 236.














Figure 1: Cyclic voltammogram of the 80% [BMIM][BF<sub>4</sub>]/H<sub>2</sub>O system. Recorded at a Pt cathode against Ag/AgCl reference. The scan rate was 0.025 Vs<sup>-1</sup>.

















































X : parts per Million : 1H

