

Initiate the camera

**Overlay the trapezoid box on the camera view
Show 3 operation buttons on the screen**

Focus button click? **Yes**

Lock auto-focus and auto-exposure

Reference button click? **Yes**

**Take a reference image; Select the areas of interest;
Save average green pixel intensities; Trigger 30 s timer**

Signal button click? **Yes**

**Take a signal image; Select the areas of interest;
Save average green pixel intensities**

**Calculate the normalized scatter intensities; Calculate and
display the bacterial concentration using standard curve**

