

Supplemental Information

for

A Paramagnetic Chemical Exchange-based MRI Probe Metabolized by Cathepsin D: Design, Synthesis and Cellular Uptake Studies

Mojmír Suchý, Robert Ta, Alex X. Li, Filip Wojciechowski, Stephen H. Pasternak, Robert Bartha and Robert H. E. Hudson*

*Correspondence to: Robert Hudson, Department of Chemistry, The University of Western Ontario, London, Ontario, Canada N6A 5B7
Fax: 1 519-661-3022; Tel: 1 519-661-2111 ext. 86349; E-mail: robert.hudson@uwo.ca

Contents:

Figure 1. Representative HPLC chromatogram of 6a	page 2
Figure 2 & 3. High resolution mass spectral characterization of 1 and 1_{cleaved}	page 3
Scheme 1. Hydrolysis of 1 by Cat D and the structure of pepstatin A.....	page 4
Figure 4. LC-MS analysis of the in-vitro metabolism of 1 and 1_{cleaved} by Cat D	page 4
Figure 5. ¹ H NMR spectrum of 4	page 5

Purification and Isolation of Compound **6a**




Figure 1: HPLC chromatogram of the crude mixture containing conjugate **6a**, $t_R = 25.1$ min.

Mass Spectral Characterization of Metalated Conjugates




Figure 2: High resolution mass spectrum (ESI-TOF) of the probe **1** (M^{5+} , left). Deconvoluted (MaxEnt 1) spectrum of **1** (right). $M=3864.6001$, conforms to the formula: $C_{166}H_{249}F_2N_{53}O_{40}STm$





Figure 3: High resolution mass spectrum (ESI-TOF) of the fragment **1_{cleaved}** (M^{3+} , left). Deconvoluted (MaxEnt 1) spectrum of **1_{cleaved}** (right). $M=1906.6001$, conforms to the formula: $C_{81}H_{101}F_2N_{15}O_{24}STm$

In-vitro Metabolism of **1** by Cathepsin D.



Scheme 1: Metabolism of the dual probe **1** by cathepsin D (Cat D); hydrolysis occurs predominantly where indicated by the asterisk. Structure of pepstatin A (**7**), an inhibitor of Cat D.




Figure 4: Total ion count spectrograms, positive mode ESI-TOF MS for the in-vitro metabolism of **1** by cathepsin D. Top spectrogram: control sample containing cathepsin D (▼) in the presence of pepstatin

A (**8**) (t_R 2.96 min). Bottom spectrogram: cleavage of probe **1** (t_R 2.19 min) to yield fragment **1_{cleaved}** (t_R 2.42 min).





Figure 5: ^1H NMR spectrum of **4** in CDCl_3 , full view (top) expanded region 1-4 ppm (bottom), showing extensive line broadening due to conformational rigidity.