

Reversible peptide particle formation using a mini amino acid sequence

Thomas B. Schuster^a, Dirk de Bruyn Ouboter^a, Enrica Bordignon^b, Gunnar Jeschke^b, Wolfgang Meier^{*a}

5

Supporting information

Materials and methods

10 Gel permeation chromatography (GPC)

Measurements of gramicidin A (gA) and Ac-X₃-gT (c = 1 mg/mL) were performed on a Shimadzu Prominence HPLC, using an Agilent PLgel column (5 μm, 10³ Å, 2 x 250 mm) and isocratic elution with THF/10% H₂O. Stabilizer-free solvent improved detection at the tryptophan adsorption maximum of 280 nm.

Tensiometer

Surface tension was measured on a Sigma 703D Tensiometer (KSV Inst., Finland) with a platinum Wilhelmy plate, pre-cleaned with ethanol and water, followed by flame annealing. Solutions were prepared 24 h prior to measurements, which were performed at room temperature.

UV-VIS

A Perkin Elmer Lambda 35 UV/VIS Spectrometer, equipped with a Peltier temperature control element and a 1 mm cuvette, was used to detect (280 nm) the onset of scattering due to particle formation.

Atomic force microscopy (AFM)

Measurements were carried out using a 5100 Agilent system (formerly PicoLE Molecular Imaging) equipped with a multi-purpose scanner. Images were acquired using a silicon cantilever (type-NCHR PointProbe® Plus, force constant 42 N/m) as indicated by the manufacturer for acoustic mode images. Samples were prepared by exposing aqueous peptide solution to a freshly cleaved mica surface for 2 min, removing the supernatant solution, and drying.

Results

Dimersiation probed by GPC

Fig. S1 illustrates GPC elution profiles of Ac-X₃-gT and gramicidin in THF and EtOH, whereby a clear separation of monomer and dimer was visible as two peaks for gramicidin only. The amphiphilic peptide elution profile was characterized by a slightly broader peak and no indication of distinct dimerisation.




Fig. S1 A) GPC elution profiles of gramicidin mono and dimer, and B) Ac-X₃-gT monomers in different solvents (THF black line, EtOH gray line)

Peptide purity and characterisation

Preparative HPLC profiles with the collected fractions for K₃-gT are presented in Fig. S2.




Fig. S2 Preparative HPLC elution profile of K₃-gT A) first run with ACN/dd H₂O 0.1% TFA B) second run using ACN/dd H₂O 2% AcOH

Verification of mass and purity of the collected fractions was analysed by MALDI-ToF-MS and HPLC. Examples are illustrated in Fig. S3.




Fig. S3 Analysis of C-K₃-gT A) HPLC profile B) MALDI-ToF-MS and AcC-X₃-gT C) HPLC profile D) MALDI-ToF-MS

60

Table S1 presents a summary of the measured mass and purity of the purified compounds.

Table S1 Verification of mass and purity of the amphiphilic peptides used.

Code	Mass g/mol	Mass (MALDI-ToF-MS)	Purity (Area %, analyt. HPLC)
K ₃ -gT	1485,9	1486.7	> 99 %
C-K ₃ -gT	1589,9	1588.6	97 %
Ac-X ₃ -gT	1653,9	1677.3	98 %
AcC-X ₃ -gT	1757,9	1779.0	96 %
AcC(sl)-X ₃ -gT	1947,0	1948.5	97 %

5

Behaviour of C-K₃-gT as surfactant

Amphiphilic peptides behave like surfactants. It therefore was possible to determine a critical micelle concentration (cmc), as illustrated in Fig. S4).




Fig. S4 Cmc of C-K₃-gT.

Self-assembly of C-K₃-gT into micelles




Fig. S5 Representative TEM image of C-K₃-gT micelles. The minor peptide bead fraction, which is not present in this micrograph, can be seen in Fig. 2 (blank)

Peptide bead formation due to charge shielding




Fig. S6 A) SEM of C-K₃-gT with NaH₂PO₄ (0,025mol/L) B) TEM of AcC-K₃-gT with NaSCN (0.1mol/L)

20

Self-assembly of Ac-X₃-gT

Self-assembled spherical structures - so-called peptide beads - from Ac-X₃-gT were imaged using SEM, TEM, and AFM (Fig. S7).




Fig. S7 A) SEM of Ac-X₃-gT peptide beads, B) TEM stained with uranyl acetate and C) AFM.

25

Cw EPR measurements in different ethanol/water mixtures




Fig. S8 cw EPR spectra of AcC-X₃-gT / AcC(sl)-X₃-gT mixtures depending of ethanol concentration. Arrows mark the characteristic features of the mobile species, which are getting dominant with increasing ethanol concentration.

10

15

Correlation functions




Fig. S9 Correlation functions of AcC-X₃-gT/AcC(sI)-X₃-gT mixture in A) 40 wt % ethanol and B) 4 wt % ethanol

5

Equilibrium of Ac-X₃-gT in ethanol mixtures

Equilibrium between micelles and peptide beads arises in ethanol/water mixtures. Converting peptide solutions at 4 wt % ethanol to 40 wt % ethanol and vice versa (using dialysis) always produces the same characteristic DLS signals, proving the reversibility of the system. A diagram of the experiment and the corresponding results are presented in Fig. S7.




Fig. S10 Reversibility of the system: Two different starting points (4 wt % and 40 wt %) yield the same DLS features at 4 wt % (■) and at 40 wt % (○) final ethanol content.

Temperature-dependent self-assembly of Ac-X₃-gT

The self-assembly process was initiated by solvent exchange from 100 w % to 20 wt % ethanol. A reversible, temperature-dependent peptide bead formation was monitored at that ethanol concentration. Temperature-dependent UV absorption measurements tracked bead formation as absorption increased due to scattering. As Fig. S11 shows, peptid beads begin to assemble below 25 °C at 20 wt % ethanol. This process was reversible. The cooling (black line) - heating (red line) cycle exhibits a slight hysteresis due to the experimental setup.

25




Fig. S11 UV measurement of temperature-dependent bead formation (heating (red line), cooling (black line)) using Ac-X₃-gT at 20 wt % EtOH. Photographs were taken at 4 °C (left) and RT (22 °C, right).

Size dependency of peptide beads in dd H₂O upon initial peptide concentration




Fig. S12 SEM micrographs of Ac-X₃-gT at A) 1.0 mol/L B) 1.5 mol/L

Notes

^a Department of Chemistry, University of Basel, Klingelbergstrasse 80, CH- 4056 Basel, Switzerland; E-mail: wolfgang.meier@unibas.ch

^b Laboratory for Physical Chemistry, ETH Zürich, Wolfgang-Pauli-Str. 10, CH- 8093 Zürich, Switzerland