Electronic Supplementary Material (ESI) for Journal of Materials Chemistry B. This journal is © The Royal Society of Chemistry 2014

## **Supporting information**

## Iridium oxide nanoparticles induced dual catalytic/inhibition based detection of phenol and pesticide compounds

Carmen C. Mayorga-Martinez<sup>a‡</sup>, Flavio Pino<sup>a‡</sup>, Sevinc Kurbanoglu<sup>a, b‡</sup>, Lourdes Rivas<sup>a</sup>, Sibel A. Ozkan<sup>b</sup>, Arben Merkoçi<sup>a,c</sup>\*

<sup>a</sup>Nanobioelectronics & Biosensors Group, ICN2- Institut Catala de Nanociencia i Nanotecnologia, Campus UAB, 08193, Bellaterra, Barcelona, Spain

<sup>b</sup>Ankara University, Faculty of Pharmacy, Department of Analytical Chemistry, 06100, Tandogan, Ankara, Turkey

<sup>c</sup>ICREA, Barcelona, Catalonia, Spain.

\* E-mail: arben.merkoci@icn.cat

‡C. C. M.-M., F. P. and S. K. contributed equally to this work.



Figure S1. Iridium 4f spectrum of IrOx NPs

Electrode	Enzyme	Linear Range (µM)	LOD (µM)	Sample	Working Potential (mV)	Ref
Gold disk electrodes (d.1.6 mm)	AChE	1.0 e <sup>-5</sup> -1.0 e <sup>-0</sup>	1.0 e <sup>-6</sup>	Spiked river water	350 (DPV)	[1]
Mini carbon paste electrode	AChE	1.0 e <sup>-5</sup> -1.0 e <sup>-0</sup>	4.0 e <sup>-6</sup>	-	50	[2]
gold electrode (d=1 mm)	AChE	5.0 e <sup>-5</sup> - 7.5 e <sup>-2</sup>	5.0 e <sup>-5</sup>	Vegetable sample	600	[3]
glassy carbon electrode	TYR	7.1 e <sup>-4</sup> - 2.85 e <sup>-2</sup>	5.7 e <sup>-4</sup>		6	[4]
SPE	TYR	1 e <sup>-2</sup> - 1 e <sup>-1</sup>	3.0 e <sup>-3</sup>	Spike tap and river water	-0.2	This work

## Table 1: Comparison between different chlorpyrifos biosensors reported in the literature.

## REFERENCES

- [1] S. Viswanathan, H. Radecka, J. Radecki, Biosens. Bioelectron. 2009, 24, 2772-2777
- [2] L. Zamifir, L. Rotariu, C. Bala, Biosens. Bioelectron. 2011, 26, 3692-3695.
- [3] C. Zhai, X. Sun, W. Sao, Z. Gong, X. Wang, Biosens. Bioelectron. 2013, 42, 124-130.
- [4] T. Liu, M. Xu, H. Yin, S. Ai, X. Qu, and S. Zong, Microchim. Acta. 2011, 175, 129-135