

Electronic Supplementary Information

Aptamer-based spectrophotometric detection of kanamycin in milk

Nandi Zhou,* Juan Zhang, Yaping Tian

The Key Laboratory of Industrial Biotechnology, Ministry of Education, School of
Biotechnology, Jiangnan University, Wuxi 214122, China

*Corresponding author. Tel: +86-510-85197831; Fax: +86-510-85197831;

E-mail: zhounandi@jiangnan.edu.cn




Fig. S1. Color change of AuNPs/ssDNA/aptamer solution. Lane 1: without any antibiotics, lane 2: with 500 nM kanamycin, lane 3: with 500 nM streptomycin, lane 4: with 500 nM tetracycline.