

Table: Details on validation of method

S. No.	Parameters	Results		
1.	R _F	0.33		
2.	Dynamic range (ng/spot)	1000-6000		
3.	Accuracy and precision	Instrumental precision	Method Precision	
		(% RSD)	(% RSD)	
		0.090	0.64	0.75
4.	% Recovery	80 % standard addition	96.24	
		100 % standard addition	97.10	
		120 % standard addition	96.60	
5.	Limit of detection (ng/spot)	12.56		
6.	Limit of quantification (ng/spot)	38.15		
7.	Linearity (correlation coefficient)	0.96105		