

Electronic supplementary information

Effects of Terminal Group and Chain Length on Temperature-Responsive Chromatography Utilizing Poly(*N*-isopropylacrylamide) Synthesized via RAFT Polymerization

Yuki Hiruta^a, Yuhei Nagumo^a, Atsushi Miki^a, Teruo Okano^b, and Hideko Kanazawa^{*,a}

^aFaculty of Pharmacy, Keio University,
1-5-30 Shibakoen, Minato, Tokyo 105-8512, Japan

^bInstitute of Advanced Biomedical Engineering and Science, Tokyo Women's Medical University,
TWIns, 8-1 Kawadacho, Shinjuku, Tokyo 162-8666, Japan

To whom correspondence should be addressed.

Email: kanazawa-hd@pha.keio.ac.jp

Table of Contents

Fig. S1	S2
Fig. S2	S2
¹ H NMR spectrum of 20000-D	S3
¹ H NMR spectrum of 20000-M	S4
¹ H NMR spectrum of 5000-D	S5
¹ H NMR spectrum of 5000-M	S6




Fig. S1 Microcalorimetric endotherm of polymers in water.






Fig. S2 Temperature dependent retention factors of steroids on HPLC column: (a) hydrocortisone, (b) prednisolone, (c) dexamethasone, (d) hydrocortisone acetate, and (e) testosterone.


¹H NMR spectrum of 20000-D



¹H NMR spectrum of 20000-M



¹H NMR spectrum of 5000-D



¹H NMR spectrum of 5000-M

