Electronic supplementary information

Illuminating trehalose mediated inhibition of protein aggregation through lysozyme-silver nanoparticle interaction

Soumik Siddhanta, a,b# Ishan Barmanb,c* and Chandrabhas Narayana*

- a) Light Scattering Laboratory, Chemistry & Physics of Materials Unit, Jawaharlal Nehru Centre for Advanced Scientific Research, Jakkur P.O., Bangalore 560 064, India .E-mail: cbhas@jncasr.ac.in
- b) Department of Mechanical Engineering, Johns Hopkins University, Baltimore, Maryland 21218, USA. E-mail: ibarman@jhu.edu
- c) Department of Oncology, Johns Hopkins University School of Medicine, Baltimore, Maryland 21287, USA.
 - # Author's present address
 - * Corresponding authors




Figure S1 (A) UV-vis spectra of Ag nanoparticles in presence of different concentrations of lysozyme. (B) The intensity of longitudinal plasmon band at 800 cm⁻¹ at different concentration points (1-4) of lysozyme which are 0, 10⁻⁴, 10⁻⁶ and 10⁻⁹ M respectively.




Figure S2 UV-vis spectra of Ag nanoparticles in presence of 1 M NaCl with intervals of 5 mins between successive spectra. Inset shows the intensity variation of 450 nm and 800 nm peaks.