

Sample	Co per cage	EY per cage
P22 _{0.25xCo/0.25xEY}	564 ± 15	583 ± 9
P22 _{0.25xCo/0.50xEY}	633 ± 4	1112 ± 17
P22 _{0.25xCo/0.75xEY}	615 ± 10	1521 ± 16
P22 _{0.25xEY/0.25xCo}	537 ± 2	775 ± 9
P22 _{0.25xEY/0.50xCo}	1067 ± 9	821 ± 3
P22 _{0.25xEY/0.75xCo}	1506 ± 21	826 ± 6
P22 _{1xCo/1xEY}	586 ± 18	1422 ± 28
P22 _{5xCo/1xEY}	3706 ± 187	893 ± 10