

In the current study

At 15 min,	before emptying		
Pepsin used	P7000, sigma -aldrich		
protein content in the HGS (g)	200g*3%	6	g
SGF added at 15 min (volume, mL)	2.5ml/min*15 min	37.5	ml
pepsin concentration in SGF	3.2g/L SGF		
pepsin content in the HGS at 15 min	pepsin concentration*SGF adde	0.12	g
protein/enzyme ratio	protein content/enzyme content	50	50 to 1

In the Gailler et al 2012' study

protein/enzyme ratio	25 to 1	The authors mentioned in the section 3.1.1.
----------------------	---------	---