Long-term delivery of alendronate through injectable tetra-PEG hydrogel to promote osteoporosis therapy

Dawei Li,‡^a Jin Zhou,‡^a Mingming Zhang,‡^b Yuanzheng Ma,^a Yanyu Yang,*^c Xue Han*^a and Xing Wang,*^{d,e}

^a The 8th Medical Center of Chinese PLA General Hospital, Beijing 100091, China.

^b PLA Strategic Support Force Characteristic Medical Center, Beijing 100101, China

^c College of Materials Science and Engineering, Zhengzhou University, Zhengzhou 450052, China


^d Beijing National Laboratory for Molecular Sciences, Institute of Chemistry, Chinese Academy of

Sciences, Beijing 100190, China

^e University of Chinese Academy of Sciences, Beijing 100049, China

[‡]These authors contributed equally to this work.

E-mail: yyyang@zzu.edu.cn; beryl282004@sina.com; wangxing@iccas.ac.cn






Fig. S2 ¹H NMR spectrum of the PEG-NHS in CDCl₃.




Fig. S3 Cell viability analysis after cultivation for 24 and 72 h.




Fig. S4 Three-dimensional fluorescent imaging after injection of tetra-PEG hydrogel.



Fig. S5 Surgery photos of rabbits in the ovariectomized group.