Electronic Supplementary Material (ESI) for Biomaterials Science. This journal is © The Royal Society of Chemistry 2020

Supplementary Information


Injectable hydrogels based on MPEG-PCL-RGD and BMSCs for bone tissue engineering

Hyun Joo Kim^{a,b}, Su Jung You^b, Dae Hyeok Yang^b, Jin Eun^c, Hae Kwan Park^c, Moon Suk Kim^d and Heung Jae Chun*,a,b,e

^aDepartment of Biomedicine & Health Sciences, ^bInstitute of Cell and Tissue Engineering, ^eDepartment of Medical Life Sciences College of Medicine, The Catholic University of Korea, Seoul 06591, Republic of Korea

^cDepartment of neurosurgery, Eunpyeong St. Mary's Hospital, College of Medicine, The Catholic University of Korea, Seoul 03312, Republic of Korea

^dDepartment of Molecular Science and Technology, Ajou University, Suwon, 16499, Republic of Korea



Supplementary Fig. 1 Evaluation of MSCs proliferation in hydrogels. Quantification of dsDNA in MP and MP/MP-RGD hydrogels by PicoGreen assay on days 0, 3, 6, 9 and 12.