

## Supporting Information

### Fluorescence Resonance Energy Transfer-based Screening for Protein Kinase C Ligands Using 6-Methoxynaphthalene-labeled 1,2-Diacylglycerol-lactones

Kohei Tsuji, Takahiro Ishii, Takuya Kobayakawa, Nami Ohashi, Wataru Nomura, and Hirokazu  
Tamamura\*

#### Table of contents

I. Fluorescent spectra in titration of compounds 4-9 .....	2
--	---

## I. Fluorescent spectra in titration of compounds 4-9

### Fluorescence spectra (ex340nm/em350-650nm)

