

Supplementary Information

Table 1. Sample weights measured before and after pyrolysis runs and solid mass yields. Weights in grams dry basis (db)

Holding Temp.	Walnut Shells			Almond Shells		
	Initial (g db)	Final (g db)	Yield (%w/w)	Initial (g db)	Final (g db)	Yield (%w/w)
300°C	0.12	0.07	41.7	0.12	0.07	58.3
350°C	0.12	0.05	58.3	0.12	0.05	41.7
450°C	0.12	0.04	33.3	0.12	0.04	33.3
550°C	0.12	0.03	25.0	0.12	0.03	25.0