Electronic Supplementary Material (ESI) for Catalysis Science & Technology. This journal is © The Royal Society of Chemistry 2023

Supporting Information

CuO/NaNbO₃ nanorod piezo-photocatalysts promoted tetracycline removal: a dynamic internal electric field

Hongtao Wang a, Haiyue Zhang a, Ziyang Long a, Haifeng Shi *ab

^a School of Science, Jiangnan University, Wuxi, 214122, P. R. China

^b National Laboratory of Solid State Microstructures, Nanjing University, Nanjing, 210093, P. R.

China

^{*}Corresponding author. E-mail: hfshi@jiangnan.edu.cn




Figure S1. XPS spectra of CuO, NNO, and CNNO-0.15 (survey).




Figure S2. The curves of the TC degradation rate constants under different conditions:

(a) Ultrasound, (b) Light, and (c) Light + Ultrasound.




Figure S3. TC degradation curve without photocatalyst under ultrasound and dark (Green), and TC degradation curve without photocatalyst under light irradiation (Red).