

Patient-convenient Long-term Alopecia Treatment via PLGA Microsphere-Loaded Candlelit Microneedle

Jiwoo Shin^a, Sung min Cho^a, Youseong Kim^a, Geonwoo Kang^b, Tobias Braun^c, Hermann Tenor^c, Christian Ludin^c, Reto Naef^{c, †}, Hyungil Jung^{a, b, †}

a. Department of Biotechnology, Yonsei University, 50 Yonsei-ro, Seodaemun-gu, Seoul, 03722, Republic of Korea

b. JUVIC Inc., 272 Digital-ro, Guro-gu, Seoul, 08389, Republic of Korea

c. TOPADUR Pharma AG, Grabenstrasse 11A 8952 Schlieren, Switzerland

† Correspondence: hijung@yonsei.ac.kr, Co-correspondence: reto.naef@topadur.com.




Figure S1. Image of all mice utilized in *in vivo* experiment from Day 1 to 28.




Figure S2. Cryotome result of pig cadaver skin after MP-CMN insertion. (scale bar: 1000 μm)




Figure S3. Image of MP-CMN residue after insertion for 15 minutes and heatmap indicating the height of residue. Color of the heatmap was based on the scale bar on top right. (scale bar: 1000 μm)