Supplementary information (ESI) for Integrative Biology

Combinatorial cell-3D biomaterials cytocompatibility screening for tissue engineering using bioinspired superhydrophobic substrates

Christiane L. Salgado, ab Mariana B. Oliveira and João F. Mano B. Oliveira and B. Oliveira an

The numerical values obtained from image-based pixel quantification and dsDNA quantification in the hydrogels for L929 cells and MC3T3-E1 cells can be observed in Figure S1 and Figure S2, respectively. The plots evidence the corresponding tendencies for image-based analysis and conventional analysis of the cell number in the hydrogels produced in the array-based platform.

Figure S1:




Figure S1 - Data for the dsDNA and DAPI quantifications for the L929 cell line in the array-based hydrogels

^a3B's Research Group – Biomaterials, Biodegradables and Biomimetics, University of Minho, Headquarters of the European Institute of Excellence on Tissue Engineering and Regenerative Medicine, AvePark, 4806-909 Taipas, Guimarães, Portugal. ^bICVS/3B's - PT Government Associate Laboratory, Braga/Guimarães, Portugal.

^{*}ICVS/3B's - PT Government Associate Laboratory, Braga/Gi *Corresponding author email: jmano@dep.uminho.pt

[‡]C.L.S. and M.B.O. contributed equally to this work.

Figure S2:




Figure S2- Data for the dsDNA and DAPI quantifications for the MC3T3-E1 cell line in the array-based hydrogels.