

Supporting Information

Magnetic Drug Carrier with Smart pH-Responsive Polymer Network Shell for Controlled Delivery of Doxorubicin

Wan-Fu Ma,^a Ke-Yi Wu,^a Jing Tang,^b Dian Li,^a Chuan Wei,^a Jia Guo,^a Shi-Long Wang^b
and Chang-Chun Wang,^{a,*}

^aState Key Laboratory of Molecular Engineering of Polymers and Department of Macromolecular Science, Laboratory of Advanced Materials, Fudan University, Shanghai 200433, People's Republic of China.

^bShanghai Key Laboratory of Cell Signaling and Diseases and School of Life Science and Technology, Tongji University, Shanghai 200092, China

*Address correspondence to ccwang@fudan.edu.cn




Fig. S1 Photographs of MCNC/PAA (A) and MCNC/PAA-DOX (B) dispersion (2 mg/mL) in serum-containing PBS before (left) and after standing for 2 h (right) at 37 °C.