



Supplementary Figure 1 Real-Time PCR Validation of Affymetrix Gene Expression Array Results . Figure 1A shows both Affymetrix fold change and qPCR fold change results for four up-regulated genes. Figure 1B shows the same for four down-regulated genes. Quantitative PCR was performed with one randomly selected nickel-transformed clone. Results for chromium clones were also analyzed by qPCR¹. Figure 1C/D: qPCR fold-change for each Ni-transformed clone for one up-regulated gene (SATB2) and one down-regulated gene (EFNA5).

¹ H. Sun, H. A. Clancy, T. Kluz, J. Zavadil, M. Costa, Comparison of Gene Expression Profiles in Chromate Transformed BEAS-2B Cells. *PLoS One* 2011, 6.