

New Journal of Chemistry


Liposomes derivatized with tetrabranched Neurotensin peptide via click chemistry reactions

Antonella Accardo, Paola Ringhieri, Diego Tesauro and Giancarlo Morelli*

Department of Pharmacy, CIRPeB, University of Naples “Federico II”, & IBB CNR, Via Mezzocannone 16, 80134 Naples, Italy

Electronic Supplementary Information

(C18)₂-Peg9-N₃ ¹H NMR spectrum recorded at 400 MHz in CDCl₃/CD₃OD 50:50 CHCl₃ internal standard:



(C18)₂-Peg9-N₃ ¹³C NMR spectrum recorded at 100.512 MHz, in CDCl₃/CD₃OD 50:50 CDCl₃ internal standard

